



Rabbit Anti-Human CGA Polyclonal Antibody (DPABH-15632)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	Recombinant Protein, antigen sequence: APDVQDCPECTLQENPFFSQPGAPILQCMGCCFSRAYPTPLRSKKTMLVQKNVTSESTCC VAKSYNRVTVMGGFKVENHTACHCSTCYHKS
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Antigen affinity purified
Conjugate	Unconjugated
Applications	IHC
Format	Liquid
Size	100 µl
Buffer	40% glycerol and PBS (pH 7.2).
Preservative	0.02% Sodium Azide
Storage	Store at +4°C for short term storage. Long time storage is recommended at -20°C. Gently mix before use. Optimal concentrations and conditions for each application should be determined by the user.

BACKGROUND

Introduction

The four human glycoprotein hormones chorionic gonadotropin (CG), luteinizing hormone (LH), follicle stimulating hormone (FSH), and thyroid stimulating hormone (TSH) are dimers consisting of alpha and beta subunits that are associated noncovalently. The alpha subunits of these hormones are identical, however, their beta chains are unique and confer biological specificity. The protein encoded by this gene is the alpha subunit and belongs to the glycoprotein hormones alpha chain family. Two transcript variants encoding different isoforms have been found for this gene. [provided by RefSeq, Nov 2011]

Keywords

CGA; glycoprotein hormones, alpha polypeptide; HCG; LHA; FSHA; GPHa; TSHA; GPHA1; CG-ALPHA; glycoprotein hormones alpha chain; FSH-alpha; LSH-alpha; TSH-alpha; lutropin alpha chain; follitropin alpha chain; thyrotropin alpha chain; choriogonadotropin alpha chain; luteinizing hormone alpha chain; chorionic gonadotrophin subunit alpha; thyroid-stimulating hormone alpha chain; follicle-stimulating hormone alpha chain; chorionic gonadotropin, alpha polypeptide; follicle-stimulating hormone alpha subunit; anterior pituitary glycoprotein hormones common subunit alpha;

GENE INFORMATION

Entrez Gene ID

[1081](#)

UniProt ID

[P01215](#)