



Anti-PLAGL2 (aa 236-344) polyclonal antibody (DPABH-15995)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	Shows weak transcriptional activatory activity.
Immunogen	Recombinant fragment corresponding to amino acids 236-344 of Human PLAGL2
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Immunogen affinity purified
Conjugate	Unconjugated
Applications	ICC/IF, WB, IHC-P
Format	Liquid
Size	100 µl
Buffer	pH: 7.20; Constituents: 59% PBS, 40% Glycerol
Preservative	0.02% Sodium Azide
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.

GENE INFORMATION

Gene Name [PLAGL2 pleiomorphic adenoma gene-like 3 \[Homo sapiens \]](#)

Official Symbol	PLAGL2
Synonyms	PLAGL2; pleiomorphic adenoma gene-like 2; ZNF900; zinc finger protein PLAGL2; C2H2-type zinc finger protein; pleiomorphic adenoma-like protein 2;
Entrez Gene ID	5326
Protein Refseq	NP_002648.1
UniProt ID	Q9UPG8
Function	RNA polymerase II core promoter proximal region sequence-specific DNA binding transcription factor activity involved in positive regulation of transcription; metal ion binding; sequence-specific DNA binding; sequence-specific DNA binding transcription factor activity