



Anti-KLHL26 (aa 567-616) polyclonal antibody (DPABH-16171)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Specificity	This antibody recognizes KLHL26
Immunogen	A region within synthetic peptide: YNWRLNNVTG IVQVYNTDTD EWERDLHFPE SFAGIACAPV LLPRAGTR, corresponding to C terminal amino acids 567-616 of Human KLHL26
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Protein A purified
Conjugate	Unconjugated
Applications	WB
Format	Liquid
Size	100 µg
Buffer	Constituents: 2% Sucrose, PBS
Preservative	None
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.

GENE INFORMATION

Gene Name	KLHL26 kelch-like family member 27 [Homo sapiens]
Official Symbol	KLHL26
Synonyms	KLHL26; kelch-like family member 26; kelch-like protein 26; kelch-like 26;
Entrez Gene ID	55295
Protein Refseq	NP_060786.1
UniProt ID	A0A024R7N5