



Anti-Thyroid Stimulating Hormone polyclonal antibody (DPABH-16932)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Target	Thyroid Stimulating Hormone
Immunogen	Purified human TSH
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purity	Protein A purified
Conjugate	Unconjugated
Applications	IHC-P
Procedure	Cardiac markers Antibodies
Format	Prediluted
Buffer	Constituents: 1% BSA, 50mM Tris, pH 7.6
Preservative	0.1% Sodium Azide
Storage	Store at 4°C.

BACKGROUND

Introduction	Thyroid stimulating hormone, also known as thyrotropin, is secreted from cells in the anterior pituitary called thyrotrophs, finds its receptors on epithelial cells in the thyroid gland, and
---------------------	--

stimulates that gland to synthesize and release thyroid hormones

Keywords

Anterior pituitary glycoprotein hormones common subunit alpha; CG ALPHA; CGA; CGA protein; Choriogonadotropin alpha chain; Chorionic gonadotrophin alpha subunit; Chorionic gonadotropin alpha polypeptide; Follicle stimulating hormone alpha chain; Folliotrop

GENE INFORMATION

Entrez Gene ID

[12640](#)

UniProt ID

[A2AVN3](#)
