



Anti-KRT18P55 (full length) polyclonal antibody (DPABH-17265)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	Highly similar to Keratin, type I cytoskeletal 18.
Immunogen	Full length protein, corresponding to amino acids 1-259 of Human FLJ40504 (NP_775895.1)
Isotype	IgG
Source/Host	Mouse
Species Reactivity	Human
Purification	Protein G purified
Conjugate	Unconjugated
Applications	WB
Format	Liquid
Size	50 µg
Buffer	Constituents: 1X PBS, pH 7.2
Preservative	None
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.

GENE INFORMATION

Gene Name	KRT18P55 keratin 18 pseudogene 56 [Homo sapiens]
-----------	--

Official Symbol	KRT18P55
Synonyms	KRT18P55; keratin 18 pseudogene 55;
Entrez Gene ID	284085