



Anti-Follicle Stimulating Hormone polyclonal antibody (DPABH-18879)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Specificity This antibody reacts with the 35 kD beta follicle stimulating hormone which, in females, is involved in maturation of ovarian follicles and estrogen secretion, and in males stimulates secretion of testosterone. This antibody reacts with FSH producing cell

Target Follicle Stimulating Hormone

Immunogen Full length native protein (purified) (Human).

Isotype IgG

Source/Host Rabbit

Species Reactivity Human

Conjugate Unconjugated

Applications IHC-FoFr, IHC-Fr, IHC-P, WB

Format Liquid

Size 500 µl

Buffer pH: 7.40; Constituent: 0.01% PBS

Preservative 0.1% Sodium Azide

Storage Store at 4°C short term (1-2 weeks). Store at -20°C or -80°C. Avoid freeze / thaw cycle.

BACKGROUND

Introduction Follicle stimulating hormone enables ovarian folliculogenesis to the antral follicle stage and is essential for Sertoli cell proliferation and maintenance of sperm quality in the testis. Members of the pituitary glycoprotein hormone family, of which FSH i

Keywords AFX1; FKH1; FKHRL1; FKHRL1P2; Forkhead box O1; Forkhead box O3; Forkhead box O4; FOXO1; FOXO2; FOXO3; FOXO4; MLLT7;
