



Anti-SPATA31E1 (aa 913-1032) polyclonal antibody (DPABH-20121)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

| | |
|----------------------------|---|
| Antigen Description | SPATA31E1 (SPATA31 subfamily E, member 1) is a protein-coding gene. An important paralog of this gene is SPATA31A1. |
| Immunogen | antigen, corresponding to amino acids 913-1032 of Human C9orf79 (Q6ZUB1). |
| Isotype | IgG |
| Source/Host | Rabbit |
| Species Reactivity | Human |
| Purification | Immunogen affinity purified |
| Conjugate | Unconjugated |
| Applications | ICC/IF, IHC-P |
| Format | Liquid |
| Size | 100 µl |
| Buffer | pH: 7.20; Constituents: 59% PBS, 40% Glycerol |
| Preservative | 0.02% Sodium Azide |
| Storage | Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles. |

GENE INFORMATION

Gene Name [SPATA31E1 SPATA31 subfamily E, member 2 \[Homo sapiens \]](#)

| | |
|------------------------|---|
| Official Symbol | SPATA31E1 |
| Synonyms | SPATA31E1; SPATA31 subfamily E, member 1; C9orf79; FAM75E1; spermatogenesis-associated protein 31E1; XXyac-YM21GA2.5; protein FAM75E1; FAM75-like protein C9orf79; family with sequence similarity 75, member E1; |
| Entrez Gene ID | 286234 |
| Protein Refseq | NP_849150.3 |
| UniProt ID | Q6ZUB1 |