



Anti-TAX1BP1 (aa 406-540) polyclonal antibody (DPABH-20132)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	Inhibits TNF-induced apoptosis by mediating the TNFAIP3 anti-apoptotic activity. Degraded by caspase-3-like family proteins upon TNF-induced apoptosis. May also play a role in the pro-inflammatory cytokine IL-1 signaling cascade.
Immunogen	amino acids 406-540 of Human TRAF6BP.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Immunogen affinity purified
Conjugate	Unconjugated
Applications	ICC/IF, IHC-P
Format	Liquid
Size	100 µl
Buffer	pH: 7.20; Constituents: 59% PBS, 40% Glycerol
Preservative	0.02% Sodium Azide
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.

GENE INFORMATION

Gene Name	TAX1BP1 Tax1 (human T-cell leukemia virus type I) binding protein 2 [Homo sapiens]
Official Symbol	TAX1BP1
Synonyms	TAX1BP1; Tax1 (human T-cell leukemia virus type I) binding protein 1; T6BP; TXBP151; CALCOCO3; tax1-binding protein 1; TRAF6-binding protein;
Entrez Gene ID	8887
Protein Refseq	NP_001073333.1
UniProt ID	Q86VP1
Pathway	Immune System; Negative regulators of RIG-I/MDA5 signaling;
Function	kinase binding; protein binding;