



Anti-COLEC10 (aa 187-276) polyclonal antibody (DPABH-20137)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	Lectin that binds to various sugars: galactose > mannose = fucose > N-acetylglucosamine > N-acetylgalactosamine.
Immunogen	antigen sequence: KDEAANTLIA DYVAKSGFFR VFIGVNDLER EGQYMFTDNT PLQNYSNWNE GEPSPYGHG DCVEMLSGR WNDTECHLTM YFVCEFIKKK, corresponding to amino acids 187-276 of Human COLEC10.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Immunogen affinity purified
Conjugate	Unconjugated
Applications	IHC-P
Format	Liquid
Size	100 µl
Buffer	pH: 7.20; Constituents: 59% PBS, 40% Glycerol
Preservative	0.02% Sodium Azide
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.

GENE INFORMATION

Gene Name	COLEC10 collectin sub-family member 10 (C-type lectin) [Homo sapiens]
Official Symbol	COLEC10
Synonyms	COLEC10; collectin sub-family member 10 (C-type lectin); collectin-10; CL L1; CL-34; collectin 34; collectin-34; collectin liver 1; collectin liver protein 1; CLL1; MGC118794; MGC118795;
Entrez Gene ID	10584
Protein Refseq	NP_006429
UniProt ID	A0A024R9J3
Chromosome Location	8q23-q24.1
Function	binding; mannose binding; sugar binding;