



Rabbit Anti-Human THSD4 Polyclonal antibody (DPABH-20213)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	THSD4 Fusion Protein expressed in E. coli
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Antigen affinity purified
Conjugate	Unconjugated
Applications	IHC, ELISA Recommended concentration: IHC: 1:50-1:500
Molecular Weight	55 kDa
Positive Control	Transfected HEK-293 cells
Format	Liquid
Size	50 µl
Buffer	PBS with 0.02% sodium azide and 50% glycerol pH 7.3.
Preservative	0.02% sodium azide
Storage	Store at -20°C. Aliquoting is unnecessary for -20°C storage.

BACKGROUND

Keywords

THSD4; thrombospondin, type I, domain containing 4; ADAMTSL6; FVSY9334; PRO34005; ADAMTSL-6; thrombospondin type-1 domain-containing protein 4; ADAMTS-like protein 6; A disintegrin and metalloproteinase with thrombospondin motifs-like protein 6;

GENE INFORMATION

Entrez Gene ID[79875](#)**UniProt ID**[Q6ZMP0](#)