



Anti-EPG5 (aa 819-947) polyclonal antibody (DPABH-20481)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	Involved in autophagy. May play a role in the degradation step of the autophagy pathway.
Immunogen	antigen sequence, corresponding to amino acids 819-947 of Human EPG5.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Immunogen affinity purified
Conjugate	Unconjugated
Applications	ICC/IF, IHC-P
Format	Liquid
Size	100 µl
Buffer	pH: 7.20; Constituents: 59% PBS, 40% Glycerol
Preservative	0.02% Sodium Azide
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.

GENE INFORMATION

Gene Name	EPG5 ectopic P-granules autophagy protein 6 homolog (C. elegans) [Homo sapiens]
-----------	---

Official Symbol	EPG5
Synonyms	EPG5; ectopic P-granules autophagy protein 5 homolog (C. elegans); HEEW1; VICIS; hEPG5; KIAA1632; ectopic P granules protein 5 homolog;
Entrez Gene ID	57724
Protein Refseq	NP_066015.2
UniProt ID	Q9BTI0