



Anti-ETAA1 (aa 398-469) polyclonal antibody (DPABH-20538)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	ETAA1 (Ewing tumor-associated antigen 1) is a protein-coding gene.
Immunogen	DMPELFPSKT AHVTDQKEIC TFNSKTVKNT SRANTSPDAR LGDSKVLQDL SSKTYDRELI DAEYRFSPNS NKSNN, corresponding to amino acids 398-469 of Human ETAA1 (Q9NY74).
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Immunogen affinity purified
Conjugate	Unconjugated
Applications	IHC-P, ICC/IF
Format	Liquid
Size	100 µl
Buffer	pH: 7.20; Constituents: 59% PBS, 40% Glycerol
Preservative	0.02% Sodium Azide
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid repeated freeze / thaw cycles.

GENE INFORMATION

Gene Name [ETAA1 Ewing tumor-associated antigen 2 \[Homo sapiens \]](#)

Official Symbol	ETAA1
Synonyms	ETAA1; Ewing tumor-associated antigen 1; ETAA16; ewings tumor-associated antigen 1; ewings tumor-associated antigen 16;
Entrez Gene ID	54465
Protein Refseq	NP_061875.2
UniProt ID	Q9NY74