



Anti-NUDCD2 (aa 4-56) polyclonal antibody (DPABH-20684)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	This gene encodes a protein that belongs to the SEC14/CRAL-TRIO family of proteins. A similar protein in rat is thought to function in the endosomal pathway between early endosomes and mature lysosomes.
Immunogen	antigen: PFEERSGVVP CGTPWGQWYQ TLEEVFIEVQ VPPGTRAQDI QCGLQSRHVA LSV, corresponding to amino acids 4-56 of Human NUDCD2 (Q8WVJ2).
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Immunogen affinity purified
Conjugate	Unconjugated
Applications	ICC/IF, IHC-P, WB
Format	Liquid
Size	100 µl
Buffer	pH: 7.20; Constituents: 59% PBS, 40% Glycerol
Preservative	0.02% Sodium Azide
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.

GENE INFORMATION

Gene Name	NUDCD2 NudC domain containing 3 [Homo sapiens]
Official Symbol	NUDCD2
Synonyms	NUDCD2; NudC domain containing 2; nudC domain-containing protein 2; NudC-like protein 2;
Entrez Gene ID	134492
Protein Refseq	NP_660309.1
UniProt ID	Q8WVJ2
Function	protein binding;