



# Anti-TBC1D12 (aa 1-56) polyclonal antibody (DPABH-20797)

This product is for research use only and is not intended for diagnostic use.

## PRODUCT INFORMATION

<b>Antigen Description</b>	May act as a GTPase-activating protein for Rab family protein(s).
<b>Immunogen</b>	antigen sequence: MVGPEDAGAC SGRNPKLLPV PAPDPVGQDR KVIRATGGFG GGVGAVEPPE EADEEE, corresponding to N terminal amino acids 1-56 of Human TBC1D12 (O60347)
<b>Isotype</b>	IgG
<b>Source/Host</b>	Rabbit
<b>Species Reactivity</b>	Human
<b>Purification</b>	Immunogen affinity purified
<b>Conjugate</b>	Unconjugated
<b>Applications</b>	IHC-P
<b>Format</b>	Liquid
<b>Size</b>	100 µl
<b>Buffer</b>	pH: 7.20; Constituents: 59% PBS, 40% Glycerol
<b>Preservative</b>	0.02% Sodium Azide
<b>Storage</b>	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.

## GENE INFORMATION

<b>Gene Name</b>	<a href="#">TBC1D12 TBC1 domain family, member 13 [ Homo sapiens ]</a>
<b>Official Symbol</b>	TBC1D12
<b>Synonyms</b>	TBC1D12; TBC1 domain family, member 12; TBC1 domain family member 12;
<b>Entrez Gene ID</b>	<a href="#">23232</a>
<b>Protein Refseq</b>	<a href="#">NP_056003.1</a>
<b>UniProt ID</b>	<a href="#">O60347</a>
<b>Function</b>	Rab GTPase activator activity;