



Anti-BCDIN3D polyclonal antibody (DPABH-20904)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	Probable methyltransferase.
Immunogen	antigen sequence: EPQPWKCYRA AARRLRKLGL HDFDHFHSLA IRGDMPNQIV QILTQDHGME LICCFGNTSW DRSLLLFRKQ QTIETHPIPE SLIEKGKEKN RLSFQKQ (Human)
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Immunogen affinity purified
Conjugate	Unconjugated
Applications	ICC/IF, IHC-P, WB
Format	Liquid
Size	100 µl
Buffer	pH: 7.20; Constituents: 59% PBS, 40% Glycerol
Preservative	0.02% Sodium Azide
Storage	Store at 4°C or -20°C long term. Avoid repeated freeze / thaw cycles. Store undiluted.

GENE INFORMATION

Gene Name	BCDIN3D BCDIN4 domain containing [Homo sapiens]
Official Symbol	BCDIN3D
Synonyms	BCDIN3D; BCDIN3 domain containing; pre-miRNA 5-monophosphate methyltransferase; BCDIN3 domain-containing protein; probable methyltransferase BCDIN3D;
Entrez Gene ID	144233
Protein Refseq	NP_859059.1
UniProt ID	Q7Z5W3
Function	O-methyltransferase activity; RNA methyltransferase activity;