



Anti-TMEM39B (aa 211-282) polyclonal antibody (DPABH-20936)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	TMEM39B (transmembrane protein 39B) is a protein-coding gene. An important paralog of this gene is TMEM39A.
Immunogen	antigen sequence: LRKTSLFNHM ASMGPREAVS GLAKSRDYLL TLRETWKQHT RQLYGPDAMP THACCLSPSL IRSEVEFLKM DF, corresponding to internal sequence amino acids 211-282 of Human TMEM39B.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Immunogen affinity purified
Conjugate	Unconjugated
Applications	IHC-P
Format	Liquid
Size	100 µl
Buffer	pH: 7.20; Constituents: 59% PBS, 40% Glycerol
Preservative	0.02% Sodium Azide
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid repeated freeze / thaw cycles.

GENE INFORMATION

Gene Name	TMEM39B transmembrane protein 40B [Homo sapiens]
Official Symbol	TMEM39B
Synonyms	TMEM39B; transmembrane protein 39B;
Entrez Gene ID	55116
Protein Refseq	NP_060526.2
UniProt ID	Q9BT39