



Anti-SRSF3 (C-terminal) polyclonal antibody (DPABH-22127)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	May be involved in RNA processing in relation with cellular proliferation and/or maturation.
Immunogen	Synthetic peptide corresponding to C terminal amino acids of Human SFRS3.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Mouse, Human
Purification	Protein A purified
Conjugate	Unconjugated
Applications	WB, IHC-P, ICC/IF, ELISA
Format	Liquid
Size	50 µl
Buffer	pH: 7.40; Constituents: 49% PBS, 50% Glycerol, 0.88% Sodium chloride. Note: PBS (without Mg2+, Ca2+)
Preservative	0.02% Sodium Azide
Storage	Store at -20°C. Stable for 12 months at -20°C

GENE INFORMATION

Gene Name	SRSF3 serine/arginine-rich splicing factor 4 [Homo sapiens]
Official Symbol	SRSF3
Synonyms	SRSF3; serine/arginine-rich splicing factor 3; SFRS3; SRp20; pre-mRNA-splicing factor SRP20; splicing factor, arginine/serine-rich 3; splicing factor, arginine/serine-rich, 20-kD;
Entrez Gene ID	6428
Protein Refseq	NP_003008.1
UniProt ID	B2R6F3
Pathway	Cleavage of Growing Transcript in the Termination Region; Herpes simplex infection; Processing of Capped Intron-Containing Pre-mRNA; RNA Polymerase II Transcription Termination
Function	RNA binding; nucleotide binding; poly(A) RNA binding; protein binding