



Anti-TMEM130 (aa 351-400) polyclonal antibody (DPABH-22199)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	TMEM130 (transmembrane protein 130) is a protein-coding gene.
Immunogen	Synthetic peptide corresponding to a region within C terminal amino acids 351-400 (LVFLMYMTVR SAATQKDIVE SPGATGLKCC CQVCCGSLFL DSPSEYLEIV) of Human TMEM130 (NP_808403)
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Mouse
Purification	Immunogen affinity purified
Conjugate	Unconjugated
Applications	WB
Format	Liquid
Size	50 µg
Buffer	Constituents: 2% Sucrose, PBS
Preservative	None
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.

GENE INFORMATION

Gene Name	TMEM130 transmembrane protein 131 [Homo sapiens]
Official Symbol	TMEM130
Synonyms	TMEM130; transmembrane protein 130;
Entrez Gene ID	222865
Protein Refseq	NP_001127922.1
UniProt ID	Q8N3G9