



Anti-FAM115C polyclonal antibody (DPABH-22860)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	FAM115C (family with sequence similarity 115, member C) is a protein-coding gene. An important paralog of this gene is FAM115A.
Immunogen	antigen sequence: AGFPGNIILN CFGLSILPQT LKAGCFPVPT PEMRSYHFRK ALSQFQAILN HENG NLEKSC LAKLRVDGAA F corresponding to amino acids of Human FAM115C.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Immunogen affinity purified
Conjugate	Unconjugated
Applications	IHC-P
Format	Liquid
Size	100 µl
Buffer	pH: 7.20; Constituents: 59% PBS, 40% Glycerol
Preservative	0.02% Sodium Azide
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.

GENE INFORMATION

Gene Name	FAM115C family with sequence similarity 116, member C [Homo sapiens]
Official Symbol	FAM115C
Synonyms	FAM115C; family with sequence similarity 115, member C; FAM139A; protein FAM115C; family with sequence similarity 139, member A;
Entrez Gene ID	285966
Protein Refseq	NP_001123497.1
UniProt ID	A6NFQ2