



Anti-CASZ1 (internal region) polyclonal antibody (DPABH-22951)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	Probable transcription factor.
Specificity	DPABH-22951 detects endogenous CASZ1 proteins, both the CASZ5 and CASZ11 isoforms.
Immunogen	A synthetic peptide corresponding to an internal region of Human CASZ1.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Immunogen affinity purified
Conjugate	Unconjugated
Applications	WB
Format	Liquid
Size	100 µg
Buffer	pH: 7.20; Constituents: 0.42% Potassium phosphate, 0.88% Sodium chloride
Preservative	0.01% Sodium Azide
Storage	Store at -20°C or lower. Aliquot to avoid repeated freezing and thawing.

GENE INFORMATION

Gene Name	CASZ1 castor zinc finger 2 [Homo sapiens]
Official Symbol	CASZ1
Synonyms	CASZ1; castor zinc finger 1; CST; SRG; CAS11; ZNF693; dJ734G22.1; zinc finger protein castor homolog 1; survival-related; castor-related protein; zinc finger protein 693; castor homolog 1, zinc finger;
Entrez Gene ID	54897
Protein Refseq	NP_001073312.1
UniProt ID	B3KRV8
Function	DNA binding; metal ion binding;