



Anti-TBC1D15 (aa 1-150) polyclonal antibody (DPABH-23055)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	Acts as a GTPase activating protein for RAB7A. Does not act on RAB4, RAB5 or RAB6.
Immunogen	Recombinant fragment, corresponding to a region within N terminal amino acids 1-150 of Human TBC1D15 (Q8TC07).
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Immunogen affinity purified
Conjugate	Unconjugated
Applications	WB
Format	Liquid
Size	50 µl
Buffer	pH: 7.00; Constituents: 0.75% Glycine, 1.21% Tris, 20% Glycerol
Preservative	None
Storage	Store at -20°C or lower. Aliquot to avoid repeated freezing and thawing.

GENE INFORMATION

Gene Name	TBC1D15 TBC1 domain family, member 16 [Homo sapiens]
------------------	--

Official Symbol	TBC1D15
Synonyms	TBC1D15; TBC1 domain family, member 15; RAB7-GAP; TBC1 domain family member 15; GAP for RAB7; GTPase-activating protein RAB7; Tre2/Bub2/Cdc16 domain family member 15; Tre-2-budding uninhibited by benzimidazole-cell division cycle 16 domain, domain family member 15;
Entrez Gene ID	64786
Protein Refseq	NP_001139685.2
UniProt ID	Q8TC07
Function	Rab GTPase activator activity; protein binding;