



Anti-CASZ1 (N-terminal) polyclonal antibody (DPABH-23095)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	Probable transcription factor.
Specificity	DPABH-23095 is expected to cross react with multiple isoforms of this protein.
Immunogen	A synthetic peptide corresponding to region near the N terminal of Human CASZ1 (CAI22571).
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Immunogen affinity purified
Conjugate	Unconjugated
Applications	WB, ELISA, IHC-P
Format	Liquid
Size	50 µg
Buffer	Constituents: 0.15M Sodium chloride, 0.02M Potassium phosphate, pH 7.2
Preservative	0.01% Sodium Azide
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid repeated freeze / thaw cycles.

GENE INFORMATION

Gene Name	CASZ1 castor zinc finger 2 [Homo sapiens]
Official Symbol	CASZ1
Synonyms	CASZ1; castor zinc finger 1; CST; SRG; CAS11; ZNF693; dJ734G22.1; zinc finger protein castor homolog 1; survival-related; castor-related protein; zinc finger protein 693; castor homolog 1, zinc finger;
Entrez Gene ID	54897
Protein Refseq	NP_001073312.1
UniProt ID	B3KRV8
Function	DNA binding; metal ion binding;