



Anti-CCDC134 polyclonal antibody (DPABH-23626)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	CCDC134 (coiled-coil domain containing 134) is a protein-coding gene.
Immunogen	Synthetic peptide corresponding to 17 aa from near the N-terminus of Human CCDC134 (NP_079097).
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Mouse, Rat, Human
Purification	Immunogen affinity purified
Conjugate	Unconjugated
Applications	ICC/IF, IHC-P, WB
Format	Liquid
Size	100 µg
Buffer	Constituents: PBS
Preservative	0.02% Sodium Azide
Storage	Store at 4°C, stable for one year. Should not be exposed to prolonged high temperatures. Avoid freeze / thaw cycles.

GENE INFORMATION

Gene Name	CCDC134 coiled-coil domain containing 135 [Homo sapiens]
Official Symbol	CCDC134
Synonyms	CCDC134; coiled-coil domain containing 134; dJ821D11.3; coiled-coil domain-containing protein 134;
Entrez Gene ID	79879
Protein Refseq	NP_079097.1
UniProt ID	Q9H6E4