



# Rabbit Anti-Human CWC27 Polyclonal Antibody (DPABH-24815)

This product is for research use only and is not intended for diagnostic use.

## PRODUCT INFORMATION

|                           |                                                                                                                                                                                                      |
|---------------------------|------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Immunogen</b>          | Recombinant Protein, antigen sequence:<br>REDQTLALLNQFKSKLTQAI AETPENDIPETEVEDDEGWMSHVLQFEDKSRKVKDASMQD<br>SDTFEIDPRNPVKNRRREESK                                                                     |
| <b>Isotype</b>            | IgG                                                                                                                                                                                                  |
| <b>Source/Host</b>        | Rabbit                                                                                                                                                                                               |
| <b>Species Reactivity</b> | Human, Rat, Mouse                                                                                                                                                                                    |
| <b>Purification</b>       | Antigen affinity purified                                                                                                                                                                            |
| <b>Conjugate</b>          | Unconjugated                                                                                                                                                                                         |
| <b>Applications</b>       | IHC, WB, ICC-IF                                                                                                                                                                                      |
| <b>Format</b>             | Liquid                                                                                                                                                                                               |
| <b>Size</b>               | 100 µl                                                                                                                                                                                               |
| <b>Buffer</b>             | 40% glycerol and PBS (pH 7.2).                                                                                                                                                                       |
| <b>Preservative</b>       | 0.02% Sodium Azide                                                                                                                                                                                   |
| <b>Storage</b>            | Store at +4°C for short term storage. Long time storage is recommended at -20°C. Gently mix before use. Optimal concentrations and conditions for each application should be determined by the user. |

## BACKGROUND

|                     |                                                                                                                                                                                                                                                                                           |
|---------------------|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Introduction</b> | PPIases accelerate the folding of proteins.                                                                                                                                                                                                                                               |
| <b>Keywords</b>     | CWC27; CWC27 spliceosome-associated protein homolog (S. cerevisiae); NY-CO-10; SDCCAG10; peptidyl-prolyl cis-trans isomerase CWC27 homolog; PPIase CWC27; PPIase SDCCAG10; antigen NY-CO-10; peptidyl-prolyl cis-trans isomerase SDCCAG10; serologically defined colon cancer antigen 10; |

## GENE INFORMATION

|                       |                        |
|-----------------------|------------------------|
| <b>Entrez Gene ID</b> | <a href="#">10283</a>  |
| <b>UniProt ID</b>     | <a href="#">Q6UX04</a> |