



Anti-FOXJ1 (aa 50-150) polyclonal antibody (DPABH-25392)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	May play an important role in cell fate determination during lung development and in spermatogenesis.
Immunogen	Synthetic peptide conjugated to KLH derived from within residues 50 - 150 of Human HFH4.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Immunogen affinity purified
Conjugate	Unconjugated
Applications	WB
Format	Liquid
Size	200 µl
Buffer	pH: 7.40; Constituent: PBS
Preservative	0.02% Sodium Azide
Storage	Store at 4°C short term (1-2 weeks). Aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.

GENE INFORMATION

Gene Name	FOXJ1 forkhead box J2 [Homo sapiens]
Official Symbol	FOXJ1
Synonyms	FOXJ1; forkhead box J1; HFH4; HFH-4; FKHL13; forkhead box protein J1; forkhead-like 13; fork head homologue 4; forkhead-related protein FKHL13; forkhead transcription factor HFH-4; hepatocyte nuclear factor 3 forkhead homolog 4;
Entrez Gene ID	2302
Protein Refseq	NP_001445.2
UniProt ID	A0A024R8P1
Function	DNA binding; DNA binding, bending; RNA polymerase II distal enhancer sequence-specific DNA binding transcription factor activity; double-stranded DNA binding