



Anti-TROVE2 (aa 1-50) polyclonal antibody (DPABH-25683)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	RNA-binding protein that binds to several small cytoplasmic RNA molecules known as Y RNAs. May stabilize these RNAs from degradation.
Immunogen	Synthetic peptide, corresponding to a region within N terminal amino acids 1-50 (MEESVNQMGP LNEKQIANSQ DGYVWQVTDM NRLHRFLCFG SEGPTYIYE) of Human TROVE2 (NP_004591.2).
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Immunogen affinity purified
Conjugate	Unconjugated
Applications	WB, IP
Format	Liquid
Size	100 µg
Buffer	Constituent: 99% Tris citrate/phosphate. Note: pH 7 to 8
Preservative	0.09% Sodium Azide
Storage	Store at 4°C.

GENE INFORMATION

Gene Name	TROVE2 TROVE domain family, member 2 [Homo sapiens]
Official Symbol	TROVE2
Synonyms	TROVE2; TROVE domain family, member 2; Sjogren syndrome antigen A2 (60kDa, ribonucleoprotein autoantigen SS A/Ro); SSA2; 60 kDa SS-A/Ro ribonucleoprotein; 60 kDa ribonucleoprotein Ro; sjogren syndrome type A antigen; gastric cancer multi-drug resistance protein; Sjogren syndrome antigen A2 (60kD, ribonucleoprotein autoantigen SS-A/Ro); RO60; SSA2; RORNP;
Entrez Gene ID	6738
Protein Refseq	NP_001035828
UniProt ID	P10155
Chromosome Location	1q31
Pathway	Systemic lupus erythematosus;
Function	RNA binding;