



Anti-CUEDC2 (C-terminal) polyclonal antibody (DPABH-26040)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	Controls PGR and ESR1 protein levels through their targeting for ubiquitination and subsequent proteasomal degradation.
Immunogen	19 amino acid peptide from near the carboxy terminus of Human CUEDC2 (Q9H467)
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Immunogen affinity purified
Conjugate	Unconjugated
Applications	WB, ELISA, IHC-P, ICC
Format	Liquid
Size	50 µg
Buffer	Constituents: PBS
Preservative	0.02% Sodium Azide
Storage	Store at 4°C short term (1-2 weeks). Aliquot and store at -20°C long term. Avoid repeated freeze / thaw cycles.

GENE INFORMATION

Gene Name	CUEDC2 CUE domain containing 3 [Homo sapiens]
Official Symbol	CUEDC2
Synonyms	CUEDC2; CUE domain containing 2; C10orf66; bA18I14.5; CUE domain-containing protein 2;
Entrez Gene ID	79004
Protein Refseq	NP_076945.2
UniProt ID	Q9H467
Function	protein binding;