



Anti-CADM4 (aa 370-384) polyclonal antibody (DPABH-26214)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	Involved in the cell-cell adhesion. Has calcium- and magnesium-independent cell-cell adhesion activity. May have tumor-suppressor activity.
Immunogen	Synthetic peptide: C-REAFLNGSDGHRKE, corresponding to C terminal amino acids 370-384 of Human Cell adhesion molecule 4 (Q8NFZ8).
Isotype	IgG
Source/Host	Goat
Species Reactivity	Human
Purification	Immunogen affinity purified
Conjugate	Unconjugated
Applications	WB, IHC-P, ELISA
Format	Liquid
Size	50 µg
Buffer	Constituents: 0.5% BSA, Tris buffered saline, pH 7.3
Preservative	0.02% Sodium Azide
Storage	Store at 4°C. Aliquot and store at -20°C long-term.

GENE INFORMATION

Gene Name	CADM4 cell adhesion molecule 5 [Homo sapiens]
Official Symbol	CADM4
Synonyms	CADM4; cell adhesion molecule 4; NECL4; TSLL2; IGSF4C; Necl-4; synCAM4; TSLC1-like 2; nectin-like 4; TSLC1-like protein 2; nectin-like protein 4; immunoglobulin superfamily member 4C; immunoglobulin superfamily, member 4C;
Entrez Gene ID	199731
Protein Refseq	NP_660339.1
UniProt ID	Q8NFZ8