



Rabbit Anti-Human CNPPD1 Polyclonal Antibody (DPABH-26357)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	Recombinant Protein, antigen sequence: RRRLYYGWDWEADCSLEELSSPVADIAVELLQKAAPSPIRRLQKKYVAHVSREACISPCA MMLALVYIERLRHRNPDLQHVSSDLFLISMMVASKYLYDEGEEEEVFNDEWGAAGGVA VPTLNALERGFLSAMDWHLYTDPREI
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Antigen affinity purified
Conjugate	Unconjugated
Applications	IHC
Format	Liquid
Size	100 µl
Buffer	40% glycerol and PBS (pH 7.2).
Preservative	0.02% Sodium Azide
Storage	Store at +4°C for short term storage. Long time storage is recommended at -20°C. Gently mix before use. Optimal concentrations and conditions for each application should be determined by the user.

BACKGROUND

Keywords

CNPPD1; cyclin Pas1/PHO80 domain containing 1; CGI-57; C2orf24; protein CNPPD1; cyclin Pas1/PHO80 domain-containing protein 1;

GENE INFORMATION

Entrez Gene ID[27013](#)**UniProt ID**[Q9BV87](#)