



Anti-TBB8B (aa 327-355) polyclonal antibody (DPABH-27137)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	Synthetic peptide conjugated to KLH, corresponding to a region within C terminal amino acids 327-355 of Human TBB8B (A6NNZ2).
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Immunogen affinity purified
Conjugate	Unconjugated
Applications	WB
Format	Liquid
Size	100 µl
Buffer	Constituent: PBS
Preservative	0.09% Sodium Azide
Storage	Store at 4°C short term (1-2 weeks). Aliquot and store at -20°C long term. Avoid repeated freeze / thaw cycles.

BACKGROUND

Introduction	Function: Tubulin is the major constituent of microtubules. It binds two moles of GTP, one at an exchangeable site on the beta chain and one at a non-exchangeable site on the alpha-chain.
---------------------	---

Similarities: Belongs to the tubulin family.

Keywords

LOC100510573; tubulin beta-8 chain B-like;
