



Anti-TBC1D7 (aa 154-203) polyclonal antibody (DPABH-27398)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	May act as a GTPase activating protein for Rab family protein(s)
Immunogen	Synthetic peptide corresponding to a region within C terminal amino acids 154-203 of mouse TBC1D7 (UniProt ID: Q9D0K0, NP_080211).
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Mouse
Purification	Immunogen affinity purified
Conjugate	Unconjugated
Applications	WB
Format	Liquid
Size	50 µg
Buffer	Constituent: 100% PBS
Preservative	None
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.

GENE INFORMATION

Gene Name [Tbc1d7 TBC1 domain family, member 7 \[Mus musculus \]](#)

Official Symbol	TBC1D7
Synonyms	TBC1D7; TBC1 domain family, member 7; TBC1 domain family member 7; 2610009C09Rik;
Entrez Gene ID	67046
Protein Refseq	NP_001239568
UniProt ID	Q9D0K0
Function	GTPase activator activity; Rab GTPase activator activity; molecular_function;