



Anti-CCNB1IP1 (aa 4-53) polyclonal antibody (DPABH-27496)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	E3 ubiquitin-protein ligase. Modulates cyclin B levels and participates in the regulation of cell cycle progression through the G2 phase. Overexpression causes delayed entry into mitosis.
Immunogen	Synthetic peptide within Mouse CCNB1IP1 aa 4-53 (N terminal). The exact sequence is proprietary. Sequence: CEDMLLCNYRKCRKLSGYAWVTACSHIFCDQHGSGEFSRSPAICPACNS Database link: NP_001104589 Run BLAST with Run BLAST with
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Mouse
Purification	Immunogen affinity purified
Conjugate	Unconjugated
Applications	WB
Format	Liquid
Size	50 µg
Buffer	Constituents: 98% PBS, 2% Sucrose
Preservative	None
Storage	Shipped at 4°C. Store at 4°C short term (1-2 weeks). Upon delivery aliquot. Store at -20°C long term. Avoid freeze / thaw cycle.

GENE INFORMATION

Gene Name	Ccnb1ip1 cyclin B1 interacting protein 1 [Mus musculus]
Official Symbol	CCNB1IP1
Synonyms	CCNB1IP1; cyclin B1 interacting protein 1; E3 ubiquitin-protein ligase CCNB1IP1; mei4; Gm288; Hei10;
Entrez Gene ID	239083
Protein Refseq	NP_001104589
UniProt ID	D3Z3K2
Function	molecular_function;