



Anti-CREG1 (aa 151-200) polyclonal antibody (DPABH-28521)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	May contribute to the transcriptional control of cell growth and differentiation. Antagonizes transcriptional activation and cellular transformation by the adenovirus E1A protein. The transcriptional control activity of cell growth requires interaction with IGF2R.
Immunogen	Synthetic peptide, corresponding to a region within C terminal amino acids 151-200 (MMSGTVTKVN KTEEDYARDS LFVRHPEMKH WPSSHNWFFA KLKISRIVWL) of Mouse CREG1 (NP_035934).
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Mouse
Purification	Immunogen affinity purified
Conjugate	Unconjugated
Applications	WB
Format	Liquid
Size	50 µg
Buffer	Constituents: 97% PBS, 2% Sucrose
Preservative	None
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid repeated freeze / thaw cycles.

GENE INFORMATION

Gene Name	Creg1 cellular repressor of E1A-stimulated genes 1 [Mus musculus]
Official Symbol	CREG1
Synonyms	CREG1; cellular repressor of E1A-stimulated genes 1; protein CREG1; Creg; AA755314;
Entrez Gene ID	433375
Protein Refseq	NP_035934
UniProt ID	O88668
Pathway	Senescence and Autophagy;
Function	FMN binding; oxidoreductase activity; transcription factor binding;