



Anti-RTN1 (N-terminal) polyclonal antibody (DPABH-28643)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	May be involved in neuroendocrine secretion or in membrane trafficking in neuroendocrine cells.
Immunogen	A synthetic peptide from the N-terminal region of Mouse Reticulon 1 conjugated to an immunogenic carrier protein.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Mouse, Rat
Purification	Whole antiserum
Conjugate	Unconjugated
Applications	WB, IHC-P, IHC-Fr
Format	Liquid
Size	100 µl
Buffer	Constituents: Whole serum
Preservative	None
Storage	Store at 4°C short term (1-2 weeks). Add glycerol to a final volume of 50%, aliquot and store at -20°C. Avoid repeated freeze / thaw cycles.

GENE INFORMATION

Gene Name	RTN1 reticulon 2 [Mus musculus]
Official Symbol	RTN1
Synonyms	RTN1; reticulon 1; Nsp; reticulon-1; neuroendocrine-specific protein;
Entrez Gene ID	104001
Protein Refseq	NP_001007597.1
UniProt ID	Q05CD8