



Anti-Neuropeptide Y (full length) polyclonal antibody (DPABH-28702)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Specificity	This antibody cross-reacts in dot blot with Peptide YY (PYY) conjugated to BSA. It shows no cross-reactivity in dot blot with Substance P, Neurokinin A, Neurokinin B, Vasoactive Intestinal Peptide (VIP), Calcitonin Gene-related Peptide (CGRP) (rat), Calcitonin, and Somatostatin conjugated to BSA, or BSA.
Immunogen	Full length Neuropeptide Y (porcine) conjugated to KLH.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Mouse, Rat, Rabbit, Chicken, Human, Monkey
Purification	Whole antiserum
Conjugate	Unconjugated
Applications	IHC-FoFr, ICC/IF, IHC-P, Dot Blot, IHC-Fr, RIA
Format	Liquid
Size	50 µl
Buffer	Constituents: Whole Serum
Preservative	15mM Sodium Azide
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.

BACKGROUND

Introduction NPY is implicated in the control of feeding and in secretion of gonadotrophin-release hormone.

Keywords NPY; neuropeptide Y; PYY4; pro-neuropeptide Y; prepro-neuropeptide Y;

GENE INFORMATION

Entrez Gene ID [4852](#)

UniProt ID [A4D158](#)
