



Anti-ZCCHC8 (aa 246-295) polyclonal antibody (DPABH-28848)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	May be involved in pre-mRNA splicing.
Immunogen	Synthetic peptide corresponding to C terminal amino acids 246-295 of Rat ZCCHC8 (NP_001099399).
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Rat
Purification	Immunogen affinity purified
Conjugate	Unconjugated
Applications	WB
Format	Liquid
Size	50 µg
Buffer	Constituents: 97% PBS, 2% Sucrose
Preservative	None
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid repeated freeze / thaw cycles.

GENE INFORMATION

Gene Name [Zcchc8 zinc finger, CCHC domain containing 8 \[Rattus norvegicus \]](#)

Official Symbol	ZCCHC8
Synonyms	ZCCHC8; zinc finger, CCHC domain containing 8; zinc finger CCHC domain-containing protein 8;
Entrez Gene ID	288661
Protein Refseq	NP_001099399
UniProt ID	D3ZUL8
Function	metal ion binding; nucleic acid binding; zinc ion binding;