



Anti-15-d-delta12,14-PGJ2 polyclonal antibody (DPABH-29523)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Specificity	This antibody detects Prostaglandin J2 (25.3 %), Prostaglandin D2 (13.13 %), D12-Prostaglandin J2 (8.89 %) and 9a,11b-Prostaglandin F2 alpha (6.93 %).
Immunogen	15-deoxy-D12,14-Prostaglandin J2 conjugated to a protein carrier molecule
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Broad
Conjugate	Unconjugated
Applications	ELISA
Format	Liquid
Size	50 µl
Buffer	PBS with 0.1% sodium azide
Preservative	0.1% Sodium Azide
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.

BACKGROUND

Keywords 15 deoxy (delta) 12 14 Prostaglandin J2; 15 deoxy Delta 12,14 Prostaglandin J2