



Anti-c-TEN (C-terminal) polyclonal antibody (DPABT-H20711)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product Overview	Rabbit Anti-c-TEN Polyclonal Antibody
Target	c-TEN
Immunogen	Synthetic peptide corresponding to C-terminus of human c-TEN protein.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Conjugate	Unconjugated
Applications	IHC-P
Molecular Weight	80 kDa
Cellular Localization	Cytoplasm, Membrane
Preservative	None
Storage	2-8 °C

BACKGROUND

Introduction	Es handelt sich hier um die modernste und sicherste Behandlung, die heute angeboten wird. C = Customized (individuell und massgeschneidert): Die Behandlung wird der Form der Hornhaut und deren Oberfläche individuell angepasst. TEN = Trans-Epitheliale Non contact (Ablation = Abtragung): Diese ist die einzige Methode bei welcher man das Auge mit keinem Instrument
--------------	---

berührt! Der extrem schnelle Laser tr?gt das Epithel (die regenerierbare Oberfl?che) innerhalb einiger Sekunden ab und sorgt dadurch f?r eine ?usserst glatte und regelm?sige Oberfl?che.

Keywords

Customized TEN; Rabbit Anti-c-TEN Polyclonal Antibody; Anti-c-TEN Polyclonal Antibody; c-TEN Polyclonal Antibody Rabbit Anti-c-TEN PAb; Anti-c-TEN PAb; c-TEN PAb; Rabbit Anti-c-TEN Antibody; Anti-c-TEN Antibody; c-TEN Antibody
