



Anti-DPP4 polyclonal antibody [Biotin] (DPABY-473)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	Dipeptidyl peptidase IV (DPPIV/CD26) is a serine exopeptidase that releases Xaa-Pro dipeptides from the N-terminus of proteins and peptides. It plays an important role in T-cell costimulation, chemokine biology, type II diabetes and tumor biology. There are at least 15 DPPIV-related human genes. Some function as enzymes, including DPP8, DPP9, fibroblast activation protein (FAP), prolyl endopeptidase (PREP), and N-acylaminoacyl-peptide hydrolase (APEH) while others (DPP6/DPPX and DPP10/DPRP3) may have lost their protease properties but acquired novel functions. For example, DPP6 is a critical component of neuronal A-type K+channels.
Specificity	Detects human DPPIV/CD26 in ELISAs and Western blots. In sandwich immunoassays,less than 0.1% cross-reactivity with recombinant mouse DPPIV/CD26 is observed.
Immunogen	Mouse myeloma cell line NS0-derived recombinant human DPPIV/CD26 . Asp34-Pro766 Accession Number Q53TN1
Isotype	IgG
Source/Host	Goat
Species Reactivity	Human
Purification	Antigen Affinity-purified
Conjugate	Biotin
Applications	Western Blot, ELISA Detection (Matched Pair)
Format	Liquid
Size	50 μg

45-1 Ramsey Road, Shirley, NY 11967, USA

Tel: 1-631-624-4882 Fax: 1-631-938-8221

Buffer	Lyophilized from a 0.2 μm filtered solution in PBS with BSA as a carrier protein.
Preservative	None
Storage	Use a manual defrost freezer and avoid repeated freeze-thaw cycles.
	12 months from date of receipt, -20 to -70 °C as supplied.
	1 month, 2 to 8 °C under sterile conditions after reconstitution.
	6 months, -20 to -70 °C under sterile conditions after reconstitution.

GENE INFORMATION

Gene Name	DPP4 dipeptidyl-peptidase 4 [Homo sapiens (human)]
Official Symbol	DPP4
Synonyms	DPP4; dipeptidyl-peptidase 4; CD26; ADABP; ADCP2; DPPIV; TP103; dipeptidyl peptidase 4; ADCP-2; DPP IV; dipeptidylpeptidase 4; dipeptidyl peptidase IV; T-cell activation antigen CD26; adenosine deaminase complexing protein 2; dipeptidylpeptidase IV (CD26,
Entrez Gene ID	<u>1803</u>
Protein Refseq	NP 001926
UniProt ID	<u>P27487</u>
Chromosome Location	2q24.3
Pathway	Incretin synthesis, secretion, and inactivation; Metabolism of proteins; Peptide hormone metabolism; Protein digestion and absorption; Synthesis, secretion, and inactivation of Glucagon-like Peptide-1 (GLP-1); Synthesis, secretion, and inactivation of Glu
Function	dipeptidyl-peptidase activity; identical protein binding; protease binding; protein binding; protein homodimerization activity; receptor binding; serine-type endopeptidase activity; serine-type peptidase activity;

Email: info@creative-diagnostics.com

Tel: 1-631-624-4882 Fax: 1-631-938-8221