



Anti-4-Nitrophenyl AHD polyclonal antibody (DPABY-785)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	4 Nitrophenyl AHD antibody was raised in sheep using Aminohydantoin-BSA as the immunogen
Source/Host	Sheep
Species Reactivity	N/A
Conjugate	Unconjugated
Applications	ELISA Pr*
Format	Liquid
Size	100 µl
Buffer	Ig fraction supplied in liquid form in 20 mM phosphate buffer, 150 mM sodium chloride and 0.09% sodium azide, pH 7.2
Preservative	0.09% Sodium Azide
Storage	Store at 2-8 °C for short term storage or -20 °C for up to 2 years

BACKGROUND

Introduction	Nitrofurans are a group of antibacterials, including nitrofurantoin, nitrofurazone, etc., that are effective against a wide range of bacteria.
Keywords	4 Nitrophenyl AHD; 1-5-nitro-2-furyl)methylideneamino]imidazolidine-2 4-dione; Chemiofuran; Cistofuran; Furadantin; Furadantine; FuralanNitrofuradantin; Macrobid; Macrofantin;

