



Anti-Ceftiofur polyclonal antibody (DPABY-811)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	Ceftiofur antibody was raised in sheep using Ceftiofur-BSA as the immunogen
Source/Host	Sheep
Species Reactivity	N/A
Conjugate	Unconjugated
Applications	ELISA Pr*
Format	Liquid
Size	100 µl
Buffer	Ig fraction supplied in liquid form in 20 mM phosphate buffer, 150 mM sodium chloride and 0.09% sodium azide, pH 7.2
Preservative	0.09% Sodium Azide
Storage	Store at 2-8 °C for short term storage or -20 °C for up to 2 years

BACKGROUND

Introduction	Ceftiofur is a third generation cephalosporin antibiotic developed exclusively for veterinary use. It has been approved worldwide for treatment of respiratory disease in swine, cattle, sheep, goats and horses, and has also been approved for foot rot and mastitis in cattle.
Keywords	Ceftiofur; Desfuroylceftiofur;