



# Anti-Gentamycin polyclonal antibody (DPABY-868)

This product is for research use only and is not intended for diagnostic use.

## PRODUCT INFORMATION

<b>Immunogen</b>	Gentamycin antibody was raised in sheep using Gentamycin-BSA as the immunogen
<b>Source/Host</b>	Sheep
<b>Species Reactivity</b>	N/A
<b>Conjugate</b>	Unconjugated
<b>Applications</b>	ELISA Pr*
<b>Format</b>	Liquid
<b>Size</b>	100 µl
<b>Buffer</b>	Ig fraction supplied in liquid form in 20 mM phosphate buffer, 150 mM sodium chloride and 0.09% sodium azide, pH 7.2
<b>Preservative</b>	0.09% Sodium Azide
<b>Storage</b>	Store at 2-8 °C for short term storage or -20 °C for up to 2 years

## BACKGROUND

<b>Introduction</b>	Gentamicin is an aminoglycoside antibiotic, used to treat many types of bacterial infections, particularly those caused by Gram-negative organisms. However, gentamicin is not used for <i>Neisseria gonorrhoeae</i> , <i>Neisseria meningitidis</i> or <i>Legionella pneumophila</i> . Gentamicin is also ototoxic and nephrotoxic, with this toxicity remaining a major problem in clinical use.
<b>Keywords</b>	Gentamycin; Garamycin; Cidomycin; Septopal;