



Anti-Lasalocid polyclonal antibody (DPABY-887)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	Lasalocid antibody was raised in sheep using Lasalocid A-BTG as the immunogen
Source/Host	Sheep
Species Reactivity	N/A
Conjugate	Unconjugated
Applications	ELISA Pr*
Format	Liquid
Size	100 µl
Buffer	Ig fraction supplied in liquid form in 20 mM phosphate buffer, 150 mM sodium chloride and 0.09% sodium azide, pH 7.2
Preservative	0.09% Sodium Azide
Storage	Store at 2-8 °C for short term storage or -20 °C for up to 2 years

BACKGROUND

Introduction	Lasalocid is able to make neutral complexes with monovalent and divalent cations and transport them through apolar phase (including lipid bilayer membranes). It can transport also big organic cations like dopamine.
Keywords	Lasalocid; Bovatec;