



# Anti-Adenovirus Polyclonal antibody (DPATB-H81048)

This product is for research use only and is not intended for diagnostic use.

## PRODUCT INFORMATION

Product Overview	Goat Anti-Adenovirus Polyclonal Antibody
Target	Adenovirus
Immunogen	Hexon from adenovirus, type 2.
Isotype	IgG
Source/Host	Goat
Species Reactivity	Adenovirus
Purification	Protein G purified
Conjugate	Unconjugated
Applications	ELISA, ICC, IF
Format	Liquid
Size	100 µl
Buffer	PBS
Preservative	0.1% Sodium Azide
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.

## BACKGROUND

**Introduction**

Adenoviruses are DNA viruses generally widespread in nature that are frequently the cause of acute upper respiratory tract infections (i.e. common colds). Forty-seven known serotypes have been isolated since they were first discovered in 1953 with 3 types known

**Keywords**

Adenovirus; Adenoviridae; Aviadenovirus; ADENOVIRUS; ADENOVIRUS F