



Anti-Allergen Jun a 3 (aa 61-160) polyclonal antibody (DPATB-H81072)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product Overview	Mouse Anti-Allergen Jun a 3 Polyclonal Antibody
Target	Allergen Jun a 3
Immunogen	Recombinant fusion protein: LAAGTASARF WGRTGCTFDA SGKGSCQTGD CGGQLSCTVS GAVPATLA EY TQSDQDYDV SLVDGFNIPL AINPTNAQCT APACKADINA VCPSELKVDG , corresponding to internal sequence amino acids 61-160 of Allergen Jun a 3 from Ashe juniper (Juniperus ashei) LAAGT
Isotype	IgG
Source/Host	Mouse
Species Reactivity	N/A
Purification	Whole antiserum
Conjugate	Unconjugated
Applications	WB
Format	Liquid
Size	100 µl
Buffer	50% Glycerol, Whole serum
Preservative	None
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.

BACKGROUND

Introduction	Allergen Jun a 3 causes an allergic reaction in human.
Keywords	Jun a3; Pathogenesis related protein; Pollen allergen Jun a 3;