



Anti-Aphrodisin (aa 107-126) polyclonal antibody (DPATB-H81088)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product Overview	Rabbit Anti-Aphrodisin Polyclonal Antibody
Target	Aphrodisin
Immunogen	Synthetic Aphrodisin Peptide (aa107-126), polyLysine-conjugated (KIFFTNKNMDRAGQET NMIV) KIFFTNKNMDRAGQETNMI Run BLAST with Run BLAST with
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	N/A
Purification	Whole antiserum
Conjugate	Unconjugated
Applications	WB
Cellular Localization	Secreted.
Format	Liquid
Size	50 µl
Buffer	Whole serum
Preservative	None
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.

BACKGROUND

Introduction

This protein is found in nasal epithelium and it binds a wide variety of chemical odorants.

Keywords

Obp 1f; OBP; OBP_RAT; Obp1f; Odorant binding protein; Odorant binding protein I f; Odorant-binding protein;
