



Anti-ASCI 2 Beta (aa 15-26) polyclonal antibody (DPATB-H81096)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product Overview	Rabbit Anti-ASCI3Beta Polyclonal Antibody
Target	ASCI3Beta
Immunogen	Synthetic peptide: CSGPGRFMAREQ, corresponding to amino acids 15-26 of ASCI3Beta CSGPGRFMAREQ Run BLAST with Run BLAST with
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	N/A
Purification	Immunogen affinity purified
Conjugate	Unconjugated
Applications	IHC-P
Positive Control	Rat dorsal root ganglia.
Format	Liquid
Size	100 µl
Buffer	PBS
Preservative	0.05% Sodium Azide
Storage	Store at +4°C short term. Aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.

BACKGROUND

Introduction	The acid-sensing ion channel family (ASIC) comprises six discrete ASIC subunits (ASIC 1A, ASIC 1B, ASIC 2A, ASIC 2B, ASIC 3, ASIC 4) that detect tissue acidosis (i.e. a decrease in pH). This acidosis phenomenon is a response to tissue injury, pain and inf
Keywords	ASCI3Beta