



Anti-BAPTA polyclonal antibody (DPATB-H81173)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product Overview	Rabbit Anti-BAPTA Polyclonal Antibody
Target	BAPTA
Immunogen	BAPTA conjugated to KLH
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	N/A
Purification	Protein A purified
Conjugate	Unconjugated
Applications	WB, IP, ICC, IF
Format	Liquid
Size	50 µg
Buffer	PBS, pH 7.4
Preservative	None
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.

BACKGROUND

Introduction

BAPTA ((1,2-bis(2-Aminophenoxy)ethane-N,N,N',N'-tetraacetic acid) is a highly selective calcium chelator. BAPTA has spawned many of the most familiar fluorescent indicators for calcium and other ions. BAPTA and its analogs including BAPTA-AM have proven m

Keywords

BAPTA, (1,2-bis(2-Aminophenoxy)ethane-N,N,N',N'-tetraacetic acid