



# Anti-Benzodiazepine polyclonal antibody (DPATB-H81179)

This product is for research use only and is not intended for diagnostic use.

## PRODUCT INFORMATION

Product Overview	Sheep Anti-Benzodiazepine Polyclonal Antibody
Target	Benzodiazepine
Immunogen	Oxazepam BgG.
Isotype	IgG
Source/Host	Sheep
Species Reactivity	N/A
Purification	Whole antiserum
Conjugate	Unconjugated
Applications	ELISA
Format	Liquid
Size	250 µl
Preservative	None
Storage	Add glycerol to a final volume of 50%, aliquot and store at -20°C. Avoid repeated freeze / thaw cycles.

## BACKGROUND

Introduction	The term benzodiazepine refers to the portion of the structure composed of a benzene ring (A)
--------------	---

fused to a seven membered diazepine ring (B). Benzodiazepine is a central nervous system depressant, a sedative hypnotic agent. Benzodiazepine receptors are located

---

**Keywords**

Benzodiazepine

---