



Anti-Biotin (full length) polyclonal antibody (DPATB-H81209)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product Overview	Rabbit Anti-Biotin Polyclonal Antibody
Target	Biotin
Immunogen	Full length protein conjugated to KLH
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	N/A
Purification	Immunogen affinity purified
Conjugate	Unconjugated
Applications	WB, cELISA, IP, ICC, IF
Format	Liquid
Size	1 mg
Buffer	50% Glycerol, PBS
Preservative	None
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.

BACKGROUND

Introduction	Biotin is a watersoluble vitamin, generally classified as a B complex vitamin, also called vitamin B4. After the initial discovery of biotin, nearly forty years of research were required to establish it as a vitamin. Biotin is required by all organisms but c
Keywords	Biotin; Vitamin B4; BIOS H; BIOEPIDERM; BIOTINUM; BIOTIN(V-H); BIOTIN; IMMOBILIZED ON DEXTRAN; (+)-BIOTIN; FACTOR S; COENZYME R; PHOTOPROBE(R) BIOTIN; RONACARE(TM) BIOTIN PLUS; Vitamin H USP24; Vitamin H(D-biotin); Vitamin H(D-biotin) USP and 98; VITAMIN